

Product Specification

20734.460 Trichloroacetic acid

Material description

CAS Number Molecular formula Molecular mass

Material

Grade

GPR RECTAPUR 76-03-9 C2HCl3O2

163.39

Characteristics **Specifications** Min. 98 % Assay Melting point 55 --> 59 °C Heavy metals (as Pb) Max. 20 ppm Ignition residue (SO4) Max. 0.1 %

Signature

We certify that this batch conforms to the specifications listed above.

This document has been produced electronically and is valid without a signature.

Anja Vanhalle, Head of Laboratory - Haasrode VWR International bvba; Geldenaaksebaan 464; BE-3001 Leuven; Belgium

For Professional use in Laboratory or Manufacturing. Not for use as an Active Pharmaceutical Ingredient or Food or Animal Feed. Suitability and intended use of the product remains the responsibility of the user

VWR International LLC, Radnor Corporate Center, Building One, Suite 200, 100 Matsonford Road, Radnor, PA 19087, USA VWR International byba, Haasrode Research Park Zone 2020, Geldenaaksebaan 464, 3001 Leuven, Belgium